

Germgel

THICK LIQUID BLEACH THAT KILLS GERMS, CLEANS AND DEODORISES

DESCRIPTION

Germgel is a thick bleach containing a high level of active chlorine, detergent and alkali to clean, deodorize and disinfect in one step.

APPLICATIONS (Usage Areas)

Germgel is ideal for use in bathrooms, toilets, change rooms, urinals and other sanitary areas. It can be used on all surfaces that are not sensitive to chlorine. It is ideal for the removal of grout stains and odours from urinals, toilets, showers, basins, baths and tiles.

DIRECTIONS FOR USE

For general sanitizing and cleaning, dilute 20ml Germgel to 1 liter water. Allow 5 minutes contact time before rinsing. For more heavily soiled areas, up to 100ml Germgel to 1 liter water can be used.

BENEFITS

Broad spectrum efficacy
High level of germicidal activity
Neutralizes odours

PRECAUTIONS

Germgel contains active chlorine.

- Do not use on wool, silk and coloured fabrics
- Avoid direct contact with eyes and skin
- Wear impervious rubber gloves
- Do not ingest product
- Keep out of the reach of children

EMERGENCY TREATMENT

Skin contact: Flush immediately with soap and water, wash contaminated clothing before reuse. If irritation develops, seek medical advice.

Eye contact: Rinse immediately with plenty of water for 15 minutes with eyelids held wide open. Seek medical advice.

Inhalation: Move person to fresh air and seek medical advice.

Ingestion: Do not induce vomiting. Give 2 glasses of water to drink and seek immediate medical attention.

STORAGE

Germgel has a shelf life of 12 months. Once dispensed use immediately. Store in a cool, dry place.

Conditions to avoid:

Temperatures above 30C

Avoid direct sunlight

Incompatible materials:

Acids, oxidizing agents.

Registration Number:

ACT29GNR529/259737/070/733

SABS Regulatory Affairs and Consumer Protection